

PRODUCT DATASHEET

**FLUOROSHIELD™ Mounting Medium With
4,6-diamidino-2-phenylindole (DAPI)****Catalog number: AR-6501-01**

AR-6501-02

AR-6501-01 AR-6501-02

Volume Volume

Component

FLUOROSHIELD™ with DAPI Mounting Medium 20 mL

100 mL

Dropper bottle, tip, and blue cap Amber bottle

Description: Fluoroshield with DAPI is an aqueous mounting medium for preserving fluorescence of tissue and cell smears.

This unique formula prevents rapid photobleaching of FITC, Texas Red, AMCA, Cy2, Cy3, Cy5, Alexa fluoro 488, Alexa fluoro 594, Green fluorescent protein (GFP), tetramethyl rhodamine, Redox. Phycoerythrin (RP-E), Phycocyanin (PC), and

Allophycocyanin (APC). Fluorescence is retained during prolonged storage at 4°C in the dark. This medium does not contain phenylenediamine, which destroys immunofluorescence of Cy dyes, RP-E, PC and APC.

The concentration of DAPI is optimized for tissues, it may be necessary to dilute this mounting medium for cells. Dilute 1:2 or

1:3 with FluoroShield (FLS) mounting medium, which does not contain DAPI.

This mounting medium is fortified with DAPI which is a counter-stain for DNA. This product is to be used in situ hybridization techniques or other methods where fluorescence of DNA staining is required. DAPI excites at 360nm and emits at 460nm, producing a blue fluorescence. RNA is also stained with DAPI.

May encounter problems with FROZEN BRAIN OR OTHER FROZEN TISSUES WITH LOTS OF FAT

Intended Use: Immunofluorescence, confocal microscopy.**Reagent: Ready to use mounting medium**

o

Storage: 2-8 C, protects from light, DO NOT FREEZE.**Refractive Index: 1.34 - 1.40 .**

Measured Refractive Index (RI) of the mounting media in liquid form was confirmed to be between 1.34-1.40. Refractive indexes change when the water solvent evaporates, and mounting media dries on slides. We do not have the means to measure the refractive indexes of dry mounting mediums; however, we expect the numbers to go higher when dried. RI of water is

1.333, RI of glass is 1.5.

Procedure

- 1. Bring the vial to room temperature.
- 2. Rinse slide to be mounted with DISTILLED OR DEIONIZED WATER, touch the edges of slide on a paper towel to

remove excess water. Place slides on a flat surface away from light.

- 3. Turn the vial upside down and open the dropper to remove any air bubbles.
- 4. Apply 1-3 drops of mounting medium directly on top of the specimen.
- 5. Let stand at room temperature for about 3-5 minutes in the dark.
- 6. Apply cover slip, carefully avoiding air bubbles.
- 7. The specimen is ready for visualization under a microscope.
- 8. One can seal the edges of the cover slip with nail polish, any organic mounting medium or our Organo mounting

medium (cat. # AR-6504). If a coverslip is not sealed, air bubbles will appear in few days.

- 9. Method for applying Coverslip: Put 1-3 drops of FLS on the specimen. After 3-4 minutes apply coverslip carefully

avoiding air bubbles. Put Kimwipes on the top of coverslip, press gently to remove excess mounting medium. With 200 micropipette add organic mounting medium to seal the edges. Incubate at 37C for one hour in the dark to dry organic o mounting medium. For long term storage it is recommended that the slide be stored in the dark at 2-8 C.

o

Removal of Coverslip: Coverslip can be removed before sealing the edges. Soak slide in warm (37 C) water for a few minutes. Carefully and slowly move the coverslip. Soak in water for an additional few minutes to remove coverslip. Rinse slide several times with warm water to remove all mounting medium. The slide can be re-stained and remounted.

Limitation and warranty: Our warranty is limited to the actual price paid for the product. We are not liable for any property damage, personnel injury, time, effort or economic loss due to our product.

Safety: This product contains sodium azide as a preservative, appropriate care should be taken in handling. National Institute of Occupational Safety and Health has warning that sodium azide can react with lead, copper, brass or solder in the plumbing system and forms hydrazoic acid in acidic condition. Discard with copious amount of water. Avoid skin and eye contact with all laboratory products. Use appropriate laboratory gear, laboratory coat, gloves and safety glasses. Do not ingest any laboratory products. This product is not approved for administration in humans or animals.

Country of Origin: USA

FluoroShield is a trademark of ImmunoBioScience corp.

For research use only; not for use in diagnostic procedures. FOR IN VITRO LABORATORY USE ONLY

“In vitro laboratory products for research”

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES.